



SPECTRO GENESIS Performance Maintenance

The SPECTRO GENESIS is the powerful alternative for sequential ICP-OES and FAA instruments and the perfect tool for routine analytical tasks. Customers value the extraordinary stability robustness and its reliability. However, preserving these qualities throughout its entire operating life isn't possible without proper maintenance.

That the breakdown risk for an analytical instrument increases with time is known, but only you can say what a breakdown will cost you. It could have wide-reaching consequences.

Take advantage of the possibility to maintain SPECTRO GENESIS's out of the box performance, while reducing the risk of costly breakdowns, with regular performance maintenance inspections within the framework of SPECTRO's AMECARE Performance Services. Based on your needs and circumstances, an individual maintenance plan is created, reflecting your quality compliance requirements and risk management strategies.

The rest is automatic. Appointments for scheduled inspections are made and the work conducted; all required maintenance and preventative measures will be carried out by authorized, factory trained and certified, service engineers that are backed by a best-of-class worldwide factory support infrastructure. The correct working order of the instrument is documented with a certificate and a test seal affixed to the instrument. The lifetime of the instrument is extended and your quality compliance documented.

SPECTRO's Performance Maintenance is a guarantee for the long-term preservation of the value and performance of your SPECTRO GENESIS.



Age of Instrument

Experience shows that, without proper maintenance, the likelihood of a breakdown increases

exponentially with the age of an instrument.

risk and extends the life of the instrument.

Regular maintenance drastically decreases this







SPECTRO GENESIS Performance Maintenance

SPECTRO GENESIS Performance Maintenance

- I ISO-conforming maintenance to prevent avertable instrument breakdown
- I Assurance of long-term analytical performance
- I Factory approved testing seal visibly affixed to the instrument for ISO audit compliance

I Detailed inspection certificate with professional assessment of all important components

You also profit from the following additional services:

- I Priority support and response times
- I Discounts on replacement parts and service for the time between agreed upon inspections
- I System backup and all necessary software updates





Ensuring peak performance and maximum return on investment over the life of your equipment.

We will deliver this by

- I Proactively maintaining your equipment to maximize availability
- I Optimizing solutions to meet your unique applications.
- I Enhancing user skills by providing access to product and application experts.

A SPECTRO GENESIS Performance Maintenance inspection includes, but is not limited to, the following:

Initial Instrument Status Evaluation

A general examination of the instrument's components and functions.

General Maintenance

Cleaning of the OPI/SPI and primary optical components, data backup of the instrument software and measurement results, installation of any instrument software updates if required

Measurements

Main voltage, argon gas flows and voltage supply, optic temperature and pressure, analytical performance

Replacement of parts and components

Replacement of filters & seals. Relay and primary optical components.

Determination of the time for the next inspection

Available for SPECTRO GENESIS FES17 SPECTRO GENESIS FEE17

www.spectro.com

GERMANY

SPECTRO Analytical Instruments GmbH Boschstrasse 10 D-47533 Kleve Tel: +49.2821.892.2100 Fax: +49.2821.892.2202 spectro.sales@ametek.com

USA

SPECTRO Analytical Instruments Inc. 91 McKee Drive Mahwah, NJ 07430 Tel: +1.800.548.5809 +1.201.642.3000 Fax: +1.201.642.3091 spectro-usa.sales@ametek.com

CHINA

AMETEK Commercial Enterprise (Shanghai) CO., LTD. Part A1, A4 2nd Floor Building No.1 Plot Section No.526 Fute 3rd Road East Pilot Free Trade Zone 200131 Shanghai Tel.: +86.21.586.851.11 Fax: +86.21.586.609.69 spectro-china.sales@ametek.com

SPECTRO operates worldwide and is present in more than 50 countries. For SPECTRO near you please visit www.spectro.com/worldwide © 2018 by AMETEK, Inc. All rights reserved, Subject to technical modifications, Rev. 0 • F-16 • Registered trademarks of SPECTRO Analytical Instruments GmbH SPECTRO: USA (3,645,267); EU (005673694); "SPECTRO": EU (009693763); iCAL: USA (3,189,726), EU (003131919); "SPECTRO GENESIS": EU (004206165); USA (3,170,644)