




## SPECTROPORT Performance Maintenance

**The SPECTROPORT delivers many advantages of SPECTRO's portable OES flagship, SPECTROTEST, in a smaller, lighter package. It accurately analyzes elements such as carbon, sulfur, phosphorus, and boron. It's as fast as an XRF, with many analyses taking only a few seconds. And it enables effortless point-and-shoot performance, to minimize operator intervention and decision-making.**



Experience shows that, without proper maintenance, the likelihood of a breakdown increases exponentially with the age of an instrument. Regular maintenance drastically decreases this risk and extends the life of the instrument.

**The breakdown risk for an analytical instrument increases with its age. Only you know how much that breakdown will cost you. It could have wide-reaching consequences: On top of the costs to repair the instrument, the instrument being down may mean that you have increased material costs, you may have penalties for missed shipments and, ultimately, you may face decreased revenues and profitability.**

Take advantage of the possibility to maintain the SPECTROPORT's out of the box performance, while reducing the risk of costly breakdowns, with regular performance maintenance inspections within the framework of SPECTRO's AMECARE Performance Services. Based on your needs and circumstances, an individual maintenance plan is created, reflecting your quality compliance requirements and risk management strategies.

The rest is automatic. Appointments for scheduled inspections are made and the work conducted; all required maintenance and preventative measures will be carried out by authorized, factory trained and certified, service engineers that are backed by a best-of-class worldwide factory support infrastructure. The correct working order of the instrument is documented with a certificate and a test seal affixed to the instrument. The lifetime of the instrument is extended and your quality compliance documented.

SPECTRO's Performance Maintenance is a guarantee for the long-term preservation of the value and performance of your SPECTROPORT.



## SPECTROPORT

### Performance Maintenance

- *ISO-conforming maintenance to prevent avertable instrument breakdown*
- *Assurance of long-term analytical performance*
- *Factory approved testing seal visibly affixed to the instrument for ISO audit compliance*
- *Detailed inspection certificate with professional assessment of all important components*

*You also profit from the following additional services:*

- *Priority support and response times*
- *Discounts on replacement parts and service for the time between agreed upon inspections*
- *System backup and all necessary software updates*

**A SPECTROPORT Performance Maintenance inspection includes, but is not limited to, the following:**

#### **Initial Instrument Status Evaluation**

A general examination of the instrument's components and functions

#### **General Maintenance**

Cleaning in and outside instrument, data backup of the instrument software and measurement results, installation of any instrument software updates if required

#### **Measurements**

Input voltage and voltage supply, iCAL 2.0 parameters before and after maintenance work, argon flow/pressure, optic temperature, analytical performance

#### **Replacement of Parts and Components**

Obligatory replacement of filters, seals, electrodes, primary optical interface components, and spark stand insulation;  
Elective replacement of mirror, windows and battery (internal/external)

#### **Determination of the Time for the Next Inspection**

Available for  
SPECTROPORT PXC01



Ensuring peak performance and maximum return on investment over the life of your equipment

We will deliver this by

- Proactively maintaining your equipment to maximize availability
- Optimizing solutions to meet your unique applications
- Enhancing user skills by providing access to product and application experts

[www.spectro.com](http://www.spectro.com)

#### **GERMANY**

SPECTRO Analytical Instruments GmbH  
Boschstrasse 10  
D-47533 Kleve  
Tel: +49.2821.892.0  
Fax: +49.2821.892.2202  
[spectro.sales@ametek.com](mailto:spectro.sales@ametek.com)

#### **U.S.A.**

SPECTRO Analytical Instruments Inc.  
91 McKee Drive  
Mahwah, NJ 07430  
Tel: +1.800.548.5809  
+1.201.642.3000  
Fax: +1.201.642.3091  
[spectro-usa.sales@ametek.com](mailto:spectro-usa.sales@ametek.com)

#### **CHINA**

AMETEK COMMERCIAL ENTERPRISE (SHANGHAI) CO., LTD.  
Part A1, A4 2nd Floor Building No.1 Plot Section,  
No.526 Fute 3rd Road East  
Pilot Free Trade Zone  
200131 Shanghai, China  
Tel.: +86.21.58685111  
Fax: +86.21.58660969  
[spectro-china.sales@ametek.com](mailto:spectro-china.sales@ametek.com)